

Easy to maintain



Cassette LM 810001

LM-LiftOut and LM-TwistOut instruments are supplied in a convenient cassette that protects both the instrument and the handler during the maintenance cycle. The cassette keeps the instrument from puncturing the sterilization pouch, so sterility is more likely to be assured. The cassette can be color-coded.

Sharpening stone LM 818002



The blades, hand-finished and made from LM-DuraGradeMAX supersteel, will stay sharp and last a long time when properly used and maintained. Sharpen the extracting instruments occasionally using this round sharpening stone as a simple and quick way to keep them in top condition and always ready to use.

Root tip instruments



Root Tip Pick LM 769-771 XSI/ES

- for the removal of root tip
- suits well in all difficult to reach areas
- available with LM-ErgoSense (ES) and LM-ErgoMax (XSI) handles

LM-RootOut LM 812210

- designed for delicate root teasing procedures
- not recommended for more substantial root tip removal
- similar to Heidbrink root tip picks
- ergonomically designed handle for a good grip



Exclusively distributed by Bisco Canada

2571 Smith Street Richmond, BC., V6X 2J1 Canada
www.biscocanada.com | info@biscocanada.com

EN: 1.800.667.8811 | FR: 1.800.211.1200 | @biscocanada

We reserve the right to make any changes

M1064EN 0215

LM Extraction

LM

feel the
difference



Extract with confidence

LMLiftOut

LMLiftOut

Luxating instruments for atraumatic extraction

LM-LiftOut instruments allow tooth extraction to be done as atraumatically as possible. This is important in order to enable rapid healing and future implant placement. The thin tip of the instrument is introduced into the periodontal space and slowly advanced toward the apex of the root while twisting the instrument gently back and forth. LM-LiftOut is not suitable for use as an elevator.



Angulations and sizes

Straight blade



Curved blade



Extract with confidence & luxate atraumatically

Familiar LM ergonomics and tactile sensitivity also in extracting instruments.

LM Feature	Benefit
Ergonomic design	- comfortable to hold and easy to rotate
Secure LM-ErgoTouch surface	- light, non-slip grip
Well-balanced	- controlled grip
Lightweight	- less strain on the hands
Clear color-coding	- easy identification
LM-DuraGradeMAX supersteel	- strong and sharp blades
Tight seam at blade juncture	- hygienic



LM extraction instruments rated as Best Product 2014.



Scan the QR for the complete evaluation!

LMLiftOut

LMLiftOut

Luxating instruments for atraumatic extraction

LM-LiftOut instruments allow tooth extraction to be done as atraumatically as possible. This is important in order to enable rapid healing and future implant placement. The thin tip of the instrument is introduced into the periodontal space and slowly advanced toward the apex of the root while twisting the instrument gently back and forth. LM-LiftOut is not suitable for use as an elevator.



Angulations and sizes

Straight blade



Curved blade



Extract with confidence & luxate atraumatically

Familiar LM ergonomics and tactile sensitivity also in extracting instruments.

LM Feature	Benefit
Ergonomic design	- comfortable to hold and easy to rotate
Secure LM-ErgoTouch surface	- light, non-slip grip
Well-balanced	- controlled grip
Lightweight	- less strain on the hands
Clear color-coding	- easy identification
LM-DuraGradeMAX supersteel	- strong and sharp blades
Tight seam at blade juncture	- hygienic



LM extraction instruments rated as Best Product 2014.




Scan the QR for the complete evaluation!

LM SlimLift

Luxating instruments designed for extremely narrow places

The LM-SlimLift instruments are especially designed to reach places you cannot go with traditional luxating instruments. It allows tooth extraction to be done as atraumatically as possible. This is important in order to enable rapid healing. The extremely fine tip of the instrument is introduced into the very narrow periodontal space. The instrument is slowly advanced towards the apex of the root while turning the instrument gently back and forth. The LM-SlimLift is not suitable for use as an elevator.



The extremely slim design of the blade allows for the most atraumatic extractions

Thin LM-SlimLift instruments - a genuinely more atraumatic solution

MD, DDS, D.Med.Sc, Oral Surgeon
Peter Jungell
Helsinki, Finland

"I use the LM-SlimLift instruments currently for nearly all extractions. I also use my older preferred luxating instruments while the use of the traditional ones has become less frequent.

LM's thin luxating instrument can enter even particularly narrow periodontal spaces. When wide instruments are used, the luxating force is spread over a wider area to the sides and deeper towards the tip of the roots so that the fracturing risk of the fragile roots, in particular, is limited considerably. The majority of extractions can be performed with luxating instruments without any need for any other operative measures.

The design and material of the LM-SlimLift instrument make it possible to get an excellent grip on the instrument. At first touch, the handles may even feel too big but in clinical



work I have noticed they fit the hand perfectly. The instrument can be held either on the palm of your hand or held only with your fingers. It is easy to pick a suitable instrument for any clinical need from the broad size selection. The wide blades are appropriate for the extraction of the lower molars while the narrow tips are perfect for instance for the extraction of the narrow roots of the upper teeth.

I can recommend the LM-SlimLift instruments for all dentists: it makes sense to abandon the old and poorly functioning instruments altogether. Too often you hear of cases where the treatment of the root has taken close to two hours."

Angulations and sizes

Straight blade

LM-SlimLift S3
LM 812630



3 mm

LM-SlimLift S4
LM 812640



4 mm

LM-SlimLift S5
LM 812650



5 mm

Curved blade

LM-SlimLift C3
LM 812830



3 mm

LM-SlimLift C4
LM 812840



4 mm

LM-SlimLift C5
LM 812850




5 mm

LM SlimLift

Luxating instruments designed for extremely narrow places

The LM-SlimLift instruments are especially designed to reach places you cannot go with traditional luxating instruments. It allows tooth extraction to be done as atraumatically as possible. This is important in order to enable rapid healing. The extremely fine tip of the instrument is introduced into the very narrow periodontal space. The instrument is slowly advanced towards the apex of the root while turning the instrument gently back and forth. The LM-SlimLift is not suitable for use as an elevator.



The extremely slim design of the blade allows for the most atraumatic extractions

Thin LM-SlimLift instruments - a genuinely more atraumatic solution

MD, DDS, D.Med.Sc, Oral Surgeon
Peter Jungell
Helsinki, Finland

"I use the LM-SlimLift instruments currently for nearly all extractions. I also use my older preferred luxating instruments while the use of the traditional ones has become less frequent.

LM's thin luxating instrument can enter even particularly narrow periodontal spaces. When wide instruments are used, the luxating force is spread over a wider area to the sides and deeper towards the tip of the roots so that the fracturing risk of the fragile roots, in particular, is limited considerably. The majority of extractions can be performed with luxating instruments without any need for any other operative measures.

The design and material of the LM-SlimLift instrument make it possible to get an excellent grip on the instrument. At first touch, the handles may even feel too big but in clinical



work I have noticed they fit the hand perfectly. The instrument can be held either on the palm of your hand or held only with your fingers. It is easy to pick a suitable instrument for any clinical need from the broad size selection. The wide blades are appropriate for the extraction of the lower molars while the narrow tips are perfect for instance for the extraction of the narrow roots of the upper teeth.

I can recommend the LM-SlimLift instruments for all dentists: it makes sense to abandon the old and poorly functioning instruments altogether. Too often you hear of cases where the treatment of the root has taken close to two hours."

Angulations and sizes

Straight blade

LM-SlimLift S3
LM 812630



3 mm

LM-SlimLift S4
LM 812640



4 mm

LM-SlimLift S5
LM 812650



5 mm

Curved blade

LM-SlimLift C3
LM 812830



3 mm

LM-SlimLift C4
LM 812840



4 mm

LM-SlimLift C5
LM 812850



5 mm

LM TwistOut

Elevators for forceful extraction

LM-TwistOut instruments are suitable for tooth extraction in situations where strong force or torque must be applied.



“This elevator from LM-Instruments is specifically designed for torque applications”

Bite Magazine, August 2014
by Dr. Davis Lakatos,
Complete Dentistry, Kilcoy, QLD
Australia

“I came across this elevator at the Australian Dental Congress in Melbourne last year. As we’re a rural practice and probably see a few more extractions than in a city area it’s a very handy tool.

I find that luxators are often quite fine and if a little too much pressure is applied, they tend to bend or break. This elevator from LM-Instruments is specifically designed for torque applications.

I use this elevator for all extractions and have the full range of straight and angled models in all sizes. It’s a very sturdy, simple tool that makes extractions much easier to perform.

I had often wondered why there wasn’t a more sturdy luxator on the market but, as it turns out, this elevator was exactly what I was looking for.”

Angulations and sizes

Straight blade

LM-TwistOut S3
LM 814230



3 mm

LM-TwistOut S4
LM 814240



4 mm

LM-TwistOut S5
LM 814250



5 mm

Curved blade

LM-TwistOut C3
LM 814430



3 mm

LM-TwistOut C4
LM 814440



4 mm

LM-TwistOut C5
LM 814450



5 mm

LM TwistOut

Elevators for forceful extraction

LM-TwistOut instruments are suitable for tooth extraction in situations where strong force or torque must be applied.



“This elevator from LM-Instruments is specifically designed for torque applications”

Bite Magazine, August 2014
by Dr. Davis Lakatos,
Complete Dentistry, Kilcoy, QLD
Australia

“I came across this elevator at the Australian Dental Congress in Melbourne last year. As we’re a rural practice and probably see a few more extractions than in a city area it’s a very handy tool.

I find that luxators are often quite fine and if a little too much pressure is applied, they tend to bend or break. This elevator from LM-Instruments is specifically designed for torque applications.

I use this elevator for all extractions and have the full range of straight and angled models in all sizes. It’s a very sturdy, simple tool that makes extractions much easier to perform.

I had often wondered why there wasn’t a more sturdy luxator on the market but, as it turns out, this elevator was exactly what I was looking for.”

Angulations and sizes

Straight blade

LM-TwistOut S3
LM 814230



3 mm

LM-TwistOut S4
LM 814240



4 mm

LM-TwistOut S5
LM 814250



5 mm

Curved blade

LM-TwistOut C3
LM 814430



3 mm

LM-TwistOut C4
LM 814440



4 mm

LM-TwistOut C5
LM 814450



5 mm

Easy to maintain



Cassette
LM 810001

LM-LiftOut and LM-TwistOut instruments are supplied in a convenient cassette that protects both the instrument and the handler during the maintenance cycle. The cassette keeps the instrument from puncturing the sterilization pouch, so sterility is more likely to be assured. The cassette can be color-coded.

Sharpening stone
LM 818002



The blades, hand-finished and made from LM-DuraGradeMAX supersteel, will stay sharp and last a long time when properly used and maintained. Sharpen the extracting instruments occasionally using this round sharpening stone as a simple and quick way to keep them in top condition and always ready to use.

LMDental

LM-Instruments Oy
PL 88 (Norrbyn rantatie 8)
FI-21601 Parainen, Finland
Tel. +358 2 4546 400
Fax +358 2 4546 444
info@lm-dental.com
www.lm-dental.com

Contact your local distributor



Root tip instruments



Root Tip Pick
LM 769-771 XSI/ES

- for the removal of root tip
- suits well in all difficult to reach areas
- available with LM-ErgoSense (ES) and LM-ErgoMax (XSI) handles

LM-RootOut
LM 812210

- designed for delicate root teasing procedures
- not recommended for more substantial root tip removal
- similar to Heidbrink root tip picks
- ergonomically designed handle for a good grip



M1064EN 0215 We reserve the right to make any changes

LM Extraction



Extract with confidence

LM

feel the
difference