

**3M** Science.  
Applied to Life.™

## 3M™ Impregum™ Super Quick Polyether Impression Material

# Capture every detail in 2 minutes? Yes, it can.




Now also  
available in  
cartridge  
delivery!



### Quick set.

With just a two-minute setting time, Impregum™ Super Quick Polyether Impression Material lives up to its name and allows for fast impressions in cases like single-unit crowns, implants, inlays/onlays and smaller bridges up to 3 units.



### 100% polyether. Trusted precision.

Impregum Super Quick materials demonstrate excellent performance in moist situations and exhibit great flow behaviour, leading to high levels of accuracy and detailed, high-precision restorations.

The nature of polyether allows you to get the full 45 seconds of working time, avoiding inaccuracies associated with premature setting times common with VPS materials.




### Make your patients feel comfortable.






An improved fresh, minty taste together with a short time in the mouth significantly increases your patients' comfort.

## Impregum™ Super Quick Polyether Impression Material

- For monophasic and 1-step impression techniques
- Perfect for impressions for single-unit crowns, implants, inlays/onlays and smaller bridges up to 3 units
- Available in heavy, medium, light and ultra light body consistencies
- Impregum Super Quick material is scannable with the current generation of scanners for a more efficient workflow



# The fast way to polyether precision for smaller cases.

	Dispensing system	Impression technique	Single-unit crown	Small implant (± 2 units)	Inlay/onlay (± 2 units)	Small bridge (± 3 units)	Recommended max. working time (min:sec)	Intra-oral setting time at 37°C (min:sec)
3M™ Impregum™ Penta™ Super Quick Polyether Impression Material (Medium Body)		Monophase or 1-step	✓	✓	✓	✓	0:45	2:00
3M™ Impregum™ Super Quick Polyether Impression Material (Medium Body)		Monophase or 1-step	✓	✓	✓	✓	0:45	2:00
3M™ Impregum™ Penta™ Super Quick Polyether Impression Material (Heavy Body)		1-step	✓	✓	✓	✓	0:45	2:00
3M™ Impregum™ Super Quick Polyether Impression Material (Light Body)		1-step	✓	✓	✓	✓	0:45	2:00
3M™ Impregum™ Super Quick Polyether Impression Material (Ultra Light Body)		1-step	✓	✓	✓	✓	0:45	2:00

## Ordering Information



3M™ Impregum™ Super Quick Polyether Impression Material Medium Body Garant™ Refills (69417)



3M™ Impregum™ Penta™ Super Quick Polyether Impression Material Medium Body/Light Body Introductory Kit (69415)



3M™ Impregum™ Penta™ Super Quick Polyether Impression Material Heavy Body/Light Body Introductory Kit (69407)

### Item no. Product Information

#### New 3M™ Garant™ Material Refills

- 69417** 3M™ Impregum™ Super Quick Polyether Impression Material Medium Body  
4 cartridges – 50 mL each, 5 3M™ Intra-oral Syringes – Purple, 5 3M™ Garant™ Mixing Tips – Purple
- 69379** 3M™ Impregum™ Super Quick Polyether Impression Material Light Body  
4 cartridges – 50 mL each, 5 3M™ Intra-oral Syringes – Purple, 5 3M™ Garant™ Mixing Tips – Purple
- 69418** 3M™ Impregum™ Super Quick Polyether Impression Material Ultra Light Body  
4 cartridges – 50 mL each, 5 3M™ Intra-oral Syringes – Purple, 5 3M™ Garant™ Mixing Tips – Purple

#### Introductory Kits for 3M™ Pentamix™ 3 Automatic Mixing Unit

- 69414** 3M™ Impregum™ Penta™ Super Quick Polyether Impression Material Medium Body  
1 base paste MB – 300 mL, 1 catalyst MB – 60 mL, 1 3M™ Impregum™ Penta™ Super Quick MB cartridge for 3M™ Pentamix™ 3 Automatic Mixing Unit, 10 3M™ Penta™ Mixing Tips – Red, 1 3M™ Penta™ Elastomer syringe, 1 3M™ Polyether Tray Adhesive, 3 3M™ Astringent Retraction Paste capsules, 6 3M™ Dual Arch Trays (3 posterior, 3 anterior trays)
- 69415** 3M™ Impregum™ Penta™ Super Quick Polyether Impression Material Medium Body/Light Body  
1 base paste MB – 300 mL, 1 catalyst MB – 60 mL, 2 3M™ Impregum™ Super Quick LB cartridges – 50 mL each, 1 3M™ Impregum™ Penta™ Super Quick MB cartridge for 3M™ Pentamix™ 3 Automatic Mixing Unit, 10 3M™ Penta™ Mixing Tips – Red, 5 3M™ Intra-oral Syringes – Purple, 5 3M™ Garant™ Mixing Tips – Purple, 5 3M™ Garant™ Intra-oral Tips – White, 1 3M™ Polyether Tray Adhesive, 3 3M™ Astringent Retraction Paste capsules, 6 3M™ Dual Arch Trays (3 posterior, 3 anterior trays)
- 69407** 3M™ Impregum™ Penta™ Super Quick Polyether Impression Material Heavy Body/Light Body  
1 base paste HB – 300 mL, 1 catalyst HB – 60 mL, 2 3M™ Impregum™ Super Quick LB cartridges – 50 mL each, 1 3M™ Impregum™ Penta™ Super Quick HB cartridge for 3M™ Pentamix™ 3 Automatic Mixing Unit, 10 3M™ Penta™ Mixing Tips – Red, 5 3M™ Intra-oral Syringes – Purple, 5 3M™ Garant™ Mixing Tips – Purple, 5 3M™ Garant™ Intra-oral Tips – White, 1 3M™ Polyether Tray Adhesive, 3 3M™ Astringent Retraction Paste capsules, 6 3M™ Dual Arch Trays (3 posterior, 3 anterior trays)

#### 3M™ Penta™ Material Refills

- 69385** 3M™ Impregum™ Penta™ Super Quick Polyether Impression Material Medium Body  
2 base pastes – 300 mL each, 2 catalysts – 60 mL each
- 69406** 3M™ Impregum™ Penta™ Super Quick Polyether Impression Material Heavy Body  
2 base pastes – 300 mL each, 2 catalysts – 60 mL each

#### Accessories

- 69387** 3M™ Impregum™ Penta™ Super Quick Polyether Impression Material Medium Body Cartridge
- 69401** 3M™ Impregum™ Penta™ Super Quick Polyether Impression Material Heavy Body Cartridge
- 71622** 48 3M™ Dual Arch Impression Trays Posterior
- 71623** 35 3M™ Dual Arch Impression Trays Anterior
- 30600** 3M™ Polyether Tray Adhesive  
1 bottle – 17 mL



[www.3M.ca/Dental](http://www.3M.ca/Dental)

Available in Canada from

**3M Canada**  
3M Oral Care  
300 Tartan Drive  
London, ON N5V 4M9  
Canada  
1-888-363-3685

**3M Deutschland GmbH**  
Dental Products  
Carl-Schurz-Straße 1  
41453 Neuss  
Germany

3M, 3M Science. Applied to Life., Garant, Impregum, Penta and Pentamix are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada. © 2020, 3M. All rights reserved. 2003-17128 E