

	Product	Dispensing System	Viscosity Low High	Setting Version	Maximum Working Time (23°C/73°F)	Maximum Intra-Oral Syringe Time (37°C/98°F)	Intra-Oral Setting Time (37°C/98°F)
sls	3M [™] Imprint [™] 4 Penta [™] Putty		Putty	Regular Set	1:30	_	2:30
Tray Materials	3M [™] Imprint [™] 4 Penta [™] Heavy 3M [™] Imprint [™] 4 Heavy		Heavy Body	Regular Set	2:00	_	2:00
Materials Ma	3M [™] Imprint [™] 4 Penta [™] Super Quick Heavy 3M [™] Imprint [™] 4 Super Quick Heavy		Heavy Body	Fast Set	1:15	_	1:15
	3M™ Imprint™ 4 Light	7		Regular Set	_	1:00	2:00
	3M [™] Imprint [™] 4 Super Quick Light	_		Fast Set	_	0:35	1:15
Wash M	3M [™] Imprint [™] 4 Regular	7		Regular Set	_	1:00	2:00
Wa	3M™ Imprint™ 4 Super Quick Regular	7		Fast Set	_	0:35	1:15

Ordering Information

	•	
Itam #	Draduat Information	

Intro Kits 71494 3M™ Imprint™ 4 Penta™ Heavy Intro Kit

Kit includes: 1 – 300 mL base paste; 1 – 60 mL catalyst; 1 – 50 mL Light cartridge; 1 – 50 mL Regular cartridge; 1 – Penta™ cartridge; 10 – Penta™ mixing tips, red; 5 – 3M™ Intra-Oral Syringes, green; 5 – mixing tips, yellow; 5 – intra-oral tips, yellow; 1 – VPS tray adhesive

71495 3M™ Imprint™ 4 Penta™ Super Quick Heavy Intro Kit

Kit includes: 1 – 300 mL base paste; 1 – 60 mL catalyst; 1 – 50 mL Super Quick Light cartridge; 1 – 50 mL Super Quick Regular cartridge; 1 – Penta™ cartridge; 10 – Penta™ mixing tips, red; 5 – 3M™ Intra-Oral Syringes, green; 5 – mixing tips, yellow; 5 – intra-oral tips, yellow; 1 – VPS tray adhesive

71536 3M™ Imprint™ 4 Penta™ Putty Intro Kit

Kit includes: 1 – 300 mL base paste; 1 – 60 mL catalyst; 1 – 50 mL Super Quick Light cartridge; 1 – 50 mL Regular cartridge; 1 – Penta™ cartridge; 10 – Penta™ mixing tips, red; 5 – 3M™ Intra-Oral Syringes, green; 5 – mixing tips, yellow; 5 – intra-oral tips, yellow; 1 – VPS tray adhesive

71534 3M™ Imprint™ 4 Intro Kit

Kit includes: 2 – 50 mL Heavy cartridges; 1 – 50 mL Light cartridge; 1 – 50 mL Regular cartridge; 10 – mixing tips, green; 10 – 3M™ Intra-Oral Syringes, green; 1 – VPS tray adhesive

71535 3M™ Imprint™ 4 Super Quick Intro Kit

Kit includes: 2 – 50 mL Super Quick Heavy cartridges; 1 – 50 mL Super Quick Light cartridge; 1 – 50 mL Super Quick Regular cartridge; 10 – mixing tips, green; 10 – 3M™ Intra-Oral Syringes, green; 1 – VPS tray adhesive

3M™ Penta™ Refills – Tray Materials Each refill includes: 2 – 300 mL base pastes; 2 – 60 mL catalysts

71484 3M™ Imprint™ 4 Penta™ Heavy

71485 3M™ Imprint™ 4 Penta™ Super Quick Heavy

71486 3M™ Imprint™ 4 Penta™ Putty

Refills - Tray Materials Each refill includes: 4 - 50 mL cartridges; 10 - mixing tips, green

71492 3M™ Imprint™ 4 Heavy

71493 3M™ Imprint™ 4 Super Quick Heavy

Item # Product Information

Economy Packs – Tray Materials Each Economy Pack includes: 10 – 50 mL cartridges

71531 3M™ Imprint™ 4 Heavy Economy Pack

71532 3M™ Imprint™ 4 Super Quick Heavy Economy Pack

Refills - Wash Materials

Each refill includes: 4 – 50 mL cartridges; 5 – mixing tips, yellow; 5 – 3M™ Intra-Oral Syringes, green

71488 3M™ Imprint™ 4 Light

71487 3M™ Imprint™ 4 Regular

71490 3M™ Imprint™ 4 Super Quick Light

71489 3M™ Imprint™ 4 Super Quick Regular

Accessories (includes: 50 – tips or syringes)
77949 3M™ Penta™ Mixing Tips – Red

71452 3M™ Garant™ VPS Mixing Tips – Yellow

71462 3M™ Garant™ VPS Intra-Oral Mixing Tips – Yellow

71450 3M[™] Garant[™] VPS Mixing Tips - Green 71506 3M[™] Intra-Oral Syringe - Green

3M

Available in Canada from:

3M Canada 300 Tartan Drive London, ON N5V 4M9 Canada 1-888-363-3685 www.3MESPE.ca

3M Deutschland GmbH Dental Products Carl-Schurz-Straße 1 41453 Neuss, Germany



Impressioning
A great fit should never be a surprise.

Learn more about Impressioning and the rest of our Simplified Crown and Bridge Solution:

3M.ca/GetSimplified

3M, 3M Science. Applied to Life., ESPE, Garant, Imprint, Penta and Pentamix are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada. © 2017, 3M. All rights reserved. 1702-07661 E