



Ultimate strength. It's that simple.

3M™ RelyX™ Ultimate Adhesive Resin Cement



Ideal for: CAD/CAM and glass ceramic restorations

Compare the simplicity of the 3M brand adhesive cementation system to other leading dental cements:


Only two components

Industry-leading bond strength with fewer working steps

When you need extra-strength bonding without the hassle of complex adhesive systems, choose RelyX™ Ultimate Adhesive Resin Cement in combination with Scotchbond™ Universal Adhesive. You'll get versatile, reliable adhesive cementation that's fast and easy to use – all with just two components.

- Proven, industry-leading bond strength
- Fast and easy procedure with only two components
- Total etch, selective etch or self-etch procedures
- First class aesthetics with natural fluorescence





When to use

RelyX™ Ultimate Adhesive Resin Cement in combination with 3M™ Scotchbond™ Universal Adhesive is indicated for final cementation of lower strength materials.

- CAD/CAM and glass ceramic restorations
- Resin-nano ceramics (e.g. Lava™ Ultimate CAD/CAM Restorative)
- · All-ceramic, composite or metal inlays, onlays, crowns and bridges
- Two 3-unit Maryland bridges and 3-unit inlay/onlay bridges*
- · All-ceramic or composite veneers
- All-ceramic, composite, or metal restorations to implant abutments
- Posts and screws
- * Excluded for patients with bruxism or periodontitis



RelyX™ Ultimate Introductory Kit (56891)



RelyX™ Ultimate Trial Kit - A1 (56894)



3M Canada 300 Tartan Drive London, ON N5V 4M9 Canada 1-888-363-3685 3M Deutschland GmbH **Dental Products**

Carl-Schurz-Straße 1 41453 Neuss, Germany

Recommended Applications (++ Highly recommended / + Recommended)

Indication	Metal/ Metal Based	Glass Ceramics (including Li. Disilicate) (e.g. e.max®, VITA Mark II, IPS Empress 2)	Oxide Ceramics (e.g. 3M" Lava" Plus High Translucency Zirconia, Brux Ziro, Procera®)	Resin Nano Ceramics (e.g. 3M [™] Lava [™] Ultimate CAD/CAM Restorative)	Resin Composites (e.g. VITA ENAMIC®)
Inlays/ Onlays	+	++	++	++	++
Crowns/ Bridges	+	++	+		++
Endodontic Posts	+		+		+
Veneers		+		+	+
Maryland Bridges	++		++		
On Implant Abutments	+	+	+		++

Ordering Information

Ordering information					
Item No.	Product	Contents			
56891	RelyX™ Ultimate Introductory Kit (incl. 4 colour shades)	(1) 5 mL bottle Scotchbond™ Universal Adhesive; (1) 3 mL etchant syringe; (25) dispensing tips for etchant syringe; (50) mixing wells; (50) microbrushes; (4) 8.5 g syringes of cement - 1 in each shade: A3, B0.5, A1 and Translucent; (30) mixing tips - regular; (10) mixing tips - wide + endo; (10) mixing tips - wide + intra-oral; (1) Technique Guide; (1) Instructions for Use			
56892	RelyX™ Ultimate Trial Kit (Translucent)	(1) 1.5 mL bottle Scotchbond [™] Universal Adhesive, (1) 3 mL etchant syringe, (25) dispensing tips for etchant syringe, (50) microbrushes, (1) 8.5 g syringe cement - Translucent, (10) mixing tips - regular, (5) mixing tips - wide, (5) intra-oral tips, (1) Technique Guide, (1) Instructions for Use			
56894	RelyX™ Ultimate Trial Kit (A1)	e (1) 1.5 mL bottle Scotchbond [™] Universal Adhesive, (1) 3 etchant syringe, (25) dispensing tips for etchant syring (50) microbrushes, (1) 8.5 g syringe cement - shade A (10) mixing tips - regular, (5) mixing tips - wide, (5) intratips, (1) Technique Guide, (1) Instructions for Use			

RelyX [™] Ultimate Adhesive Resin Cement Refills				
Item No.	Shade	Contents		
56887	A1	(1) 8.5 g syringe cement		
56888	A3	(10) mixing tips - regular		
56889	B0.5 Bleach	(5) mixing tips - wide (5) intra-oral tips		
56890	Translucent			

Scotchbond [™] Universal Adhesive				
Item No.	Product	Size		
41258	Scotchbond™ Universal Adhesive Refill	(1) 5 mL vial		
41263	Scotchbond™ Universal Etchant – Etching Gel	(2) 3 mL syringes		



3M, 3M Science. Applied to Life., Lava, ESPE, RelyX and Scotchbond are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada. All other trade names referenced are the service marks, trademarks or registered trademarks of their respective companies. © 2017, 3M. All rights reserved. 1612-07239 E