



**ENDODONTIC  
PRODUCTS  
AVAILABLE!**

---



**ENDODONTIC BURS**

---



**RUBBER DAM**

---



**MTA CEMENTS**



## GET THE PERFECT MIX EVERY TIME

MTA UNIVERSAL CEMENT - A MATERIAL WITH A MEDIUM CONSISTENCY SUITABLE FOR ALL INDICATIONS.

Zendo MTA is an endodontic repair cement delivered in easy to use capsules.

The MTA powder consists of very fine hydrophilic particles of several mineral oxides. After contact with MTA liquid, it forms a gel that hardens to an impermeable barrier.

Zendo MTA is delivered in convenient capsules. The capsules are easily activated and the content of the capsule is easily extruded using the Zendo Capsule Dispenser.

Capsule mixing (mixing time 30 seconds) is achieved by a high-frequency mixer (e.g. amalgamator) with about 4,300 oscillations/minute.

Zendo MTA Universal was evaluated by 21 CR evaluators and 86% said they would incorporate MTA capsules into their practice. 90% rated it excellent or good and worthy of trial by colleagues.

## Zendo MTA Universal





# MTA UNIVERSAL CEMENT


## Indications:

**Repair of Root Perforations**  
• MTA Firm NORMAL SET

**Pulp Capping**  
• MTA Flow FAST SET

**Apexification (orthograde)**  
• MTA Universal

**Root End Filling (retrograde)**  
• MTA Universal



## FEATURES:

- Biocompatible and bioactive permanent cement
- Excellent marginal adaptation due to bonding to dentin
- Very low solubility
- Bacteriostatic
- High radiopacity
- Excellent handling
- Excellent dimensional stability

## BENEFITS:


- Very easy to use
- Convenient mix and fixed proportions with capsule
- Real biological sealing
- Excellent barrier against bacteria and fluids
- Easy clinical documentation with radiographs
- Excellent marginal sealing
- Excellent choice for the non-surgical treatment of furcal and radicular strip perforations

# MTA UNIVERSAL CEMENT

A material with a medium consistency suitable for all indications.


## Clinical Case

The radiographs show a tooth with a perforated root canal. The setting reaction takes place without shrinkage and thus results in an excellent physical sealing. It follows a biological sealing by building of a protective layer against bacteria with the formation of calcium hydroxide and apatite




### BEFORE TREATMENT WITH MTA

Perforated root canal and demineralized bone matrix



### AFTER MTA TREATMENT

Repaired root canal thanks to MTA cement



### 1 YEAR AFTER TREATMENT

New bone augmentation for better retention of the tooth thanks to MTA cement



## SILK RUBBER DAM

Medium thickness | Colours - Dark



Dental dams are recommended in most of the provincial guidance documents for reducing infectious aerosol production during dental procedures.



52 5x5 sheets



36 6x6 sheets

### Key Features

- Zendo dental dam features strength, tear-resistance and durability that virtually eliminates failure.
- Pre-cut latex dental dams provide a barrier against moisture and also create an uncontaminated working field free of visual obstructions.



## ENDODONTIC BURS

Zendo Select Carbide Burs are made by one of the most reputable bur manufacturers in the industry whose award winning burs are used around the world.

Safe-End Tip NOT Coated burs 220-014M	ROUND sterile diamond burs 697-021C	SAFE END sterile diamond burs 220-014M	SELECT Carbide Burs Tapered Cone Safe End Tip 152-016	SELECT Carbide Burs Tapered Cone ZKYRA 151L-Z (long)	Zendo ZKRYA 10 sterile diamond burs 151-Z	TAPERED CONE ROUNDED TIP sterile diamond burs E164-015C	UNCOATED TIP sterile diamond burs E165-014LC	TAPERED CONE ROUNDED TIP sterile diamond burs E164-011C
<b>931113015</b>	<b>931197021</b>	<b>931120014</b>	<b>931152016</b>	<b>931115112</b>	<b>931115120</b>	<b>931164015</b>	<b>931165015</b>	<b>931164011</b>
								