

Now our most popular flowable is easier to inject.

3M™ Filtek™ Supreme Ultra Flowable Restorative is available in a new ergonomic-designed syringe. Designed for improved comfort and easier injection, this new product is called 3M™ Filtek™ Supreme Flowable Restorative.



Key benefits

- Virtually no bubbles or material "run-on"
- Easy-to-bend cannula for improved access
- Ergonomic-design is easy to hold and inject
- Excellent adaptation, polish retention and wear resistance
- Versatile for a wide range of indications and techniques

Can be warmed for injection moulding.

No higher than 70°C (158°F), no longer than 1 hour; up to 25 times.



Clinical Cases

3M™ Filtek™ Supreme Flowable Restorative offers a smooth transition from material to tooth for esthetic restorations.







Flowable placement

Photos courtesy of Dr. Gunnar Reich, Munich, Germany.

3M™ Filtek™ Supreme Flowable Restorative adapts easily and seals the margin.








Flowable placement

Photos courtesy of Dr. Giuseppe Chiodera, Brescia, Italy.



Ordering Information



3M™ Filtek™ Supreme Flowable **Restorative Refills**

Syringe refills include: 2 - 2 g syringes; 20 dispensing tips.

Shade	Item No.
A1	6032A1
A2	6032A2
A3	6032A3
A3.5	6032A3.5
A4	6032A4
B1	6032B1
B2	6032B2
C2	6032C2
D2	6032D2
White	6032W
Extra White	6032XW
OA3	6032AO3

	il .
Item No.	Product Info.
3700T	3M™ Filtek™ Supreme
	Flowable Restorative
	Dispensing Tips, 20-gauge,
	Blue (20)



Available in Canada from

3M Canada 3M Oral Care 300 Tartan Drive London, ON N5V 4M9 Canada 1-888-363-3685 3M.ca/dental

3M ESPE Dental Products 2510 Conway Avenue St. Paul, MN 55144 USA 1-800-634-2249

¹ Including small, non stress-bearing occlusal restorations.

² Repair of acrylic temporary and composite resin.